

CURA[^] FLEX

OUT-PATIENT SOLUTIONS



TABLE OF CONTENTS

01.	Electric Examination Couch	P. 03
02.	2-Function Electric Examination Couch	P. 04
03.	3-Section Electric Examination Couch	P. 05
04.	Fixed Examination Couch	P. 06
05.	Manual Examination Couch	P. 06
06.	Electric Gynaecology Examination Couch	P. 07
07.	Manual Gynaecology Couch	P. 08
08.	4-Function Hemodyalisis Chair	P. 10
09.	3-Function Trendelenburg Hemodyalisis Chair	P. 11
10.	3-Function Hemodyalisis Chair	P. 12
11.	Manual 2-Function Hemodyalisis Chair	P. 13
12.	Electric 2-Function Hemodyalisis Chair	P. 13
13.	Multi-Functional Chair	P. 14
14.	Instruments' Tables	P. 15
15.	Ward Screens Bowl Stands Steps	P. 16
16.	Veterinary Examination Couch	P. 17



ELECTRIC EXAMINATION COUCH

MLC-05



Length	Width	Height	Backrest Tilting
1850 mm	600/650 mm	570-900 mm	80°

- Equipped with a Paper Roll Holder
- Height Adjustable by Electric Foot Switch
- Optional Accessories:
 - I.V. Stand
 - Headrest Opening

2-FUNCTION ELECTRIC EXAMINATION COUCH

MLC-05S



Length	Width	Height	Backrest Tilting
1850 mm	600/650 mm	570-900 mm	80°

- Equipped with a Paper Roll Holder
- Backrest and Height Adjustable by Electric Hand Control
- Optional Accessories:
 - I.V. Stand
 - Headrest Opening

3-SECTION ELECTRIC EXAMINATION COUCH

MLC-04



Foot Switch

Length	Width	Height	Backrest Tilting
1750 mm	600/650 mm	600-950 mm	80°

- Equipped with a Paper Roll Holder
- Height Adjustable by Electric Foot Switch
- Optional Accessories:
 - I.V. Stand
 - Headrest Opening

EXAMINATION COUCHES

FIXED EXAMINATION COUCH

MLC-01E



Length	Width	Height
1850 mm	600 mm	820 mm

- Easily sanitized to guarantee longer life span
- Accessories:
 - Paper Roll Holder

MANUAL EXAMINATION COUCH

MLC-01



Length	Width	Height
1850 mm	600/650 mm	700/820 mm

- Easily sanitized to guarantee longer life span
- Manual Backrest Tilting
- Accessories:
 - Paper Roll Holder
 - I.V. Pole

ELECTRIC GYNAECOLOGY EXAMINATION COUCH

MLC-06



Adjustable
Leg Holders



S.S Basin

Foot Switch



Length	Width	Height
1750 mm	600/650 mm	600-950 mm

- Adjustable back and knee rest tilting
- Adjustable height using foot switch
- Can be extended into a full bed
- Optional Accessory:
 - Side Rails
 - I.V. Pole


MANUAL GYNAECOLOGY COUCH

MLC-02



Length	Width	Height
1750 mm	600/650 mm	700/820 mm

- Easily Sanitized to Guarantee Longer Life Span
- Manual Back and Knee Rest Tilting
- Can be Extended into a Full Bed
- Optional Accessory:
 - I.V. Pole
 - Paper Roll Holder



4-FUNCTION HEMODIALYSIS CHAIR

MLH-06I

Padded
Seat Cushion

Over-Seat
Table



Length	Width	Height	Trendelenburg
1950 mm	600 mm	600-800 mm	16°



- Can be Extended into a Full Bed
- Movable Armrests
- Adjustable Footrest
- Provided with Lower Shelf
- **Optional Accessories:**
 - Paper Roll Holder
 - Over-Seat Table

3-FUNCTION TRENDLENBURG HEMODIALYSIS CHAIR

MLH-06T



Length	Width	Height	Trendelenburg
1950 mm	600 mm	550 mm	16°



- Can be Extended into a Full Bed
- Emergency Trendelenburg Function
- Movable Armrests
- Adjustable Footrest
- Provided with Lower Shelf
- **Optional Accessories:**
 - Paper Roll Holder
 - Over-Seat Table

3-FUNCTION HEMODIALYSIS CHAIR

MLH-06



Length	Width	Height
1950 mm	600 mm	550-750 mm

- Can be Extended into a Full Bed
- Movable Armrests
- Adjustable Footrest
- Adjustable Height
- Provided with Lower Shelf
- Optional Accessories:
 - Over Bed Table
 - Paper Roll Holder

HEMODIALYSIS CHAIRS

MANUAL 2-FUNCTION HD CHAIR

MLH-15



Length	Width	Height
1950 mm	600 mm	550 mm

- Provided with an Adjustable Stainless Steel I.V Pole
- Noiseless Castors with Brakes
- Optional Accessories:
 - Over Seat Table

ELECTRIC 2-FUNCTION HD CHAIR

MLH-07



Length	Width	Height
1950 mm	600 mm	550 mm

- Can be Extended into a Full Bed
- Movable Armrests
- Adjustable Footrest
- Provided with Lower Shelf
- Optional Accessories:
 - Paper Roll Holder
 - Over Seat Table

MULTI-FUNCTIONAL CHAIR

MLC-03



Length	Width	Height
1750 mm	550 mm	470 mm

- Can be Extended into a Full Bed
- Movable Stainless Steel Armrests
- Optional Accessories:
 - I.V Infusion Stand
 - Arm Holder

INSTRUMENTS' TABLES

DUAL-TIER TABLE

MLC-31



Length	Width	Height
400/500 mm	600/700 mm	800 mm

- Stainless Steel Structure
- Moves on a 50mm Diameter Castors

TRIPLE-TIER TABLE

MLC-32



Length	Width	Height
400/500 mm	600/700 mm	800 mm

- Stainless Steel Structure
- Moves on a 50mm Diameter Castors

MAYO TABLE

MLC-34



Length	Width	Height
460 mm	630 mm	830-1170 mm

- Steel Structure with Chrome Plating
- Stainless Steel Removable Tray
- Adjustable Height using a Hand Knob
- 75mm Diameter Castors with Brakes

WITH DRAWERS

MLC-31D



Length	Width	Height
400/500 mm	600/700 mm	800 mm

- Stainless Steel Structure
- 2 Stainless Steel Drawers
- Moves on a 50mm Diameter Castors

WARD SCREENS

MLC-61



Length	Width
1800 mm	1500 mm

- Chrome-Plated Steel Tubular Frame
- Available in Distiguated Fabric (MLC-61L)
- Available in 4-Section (MLC-61Q)

MLC-62

BOWL STANDS

MLC-63



- Chrome-Plated Steel
- Easily Removed and Installed
- Heavy Duty Base with Castors



MLC-64

STEPS

MLC-65

SINGLE



Leather-Covered
Anti-Slip Protection

Length	Width	Height
400 mm	400 mm	250 mm

DOUBLE



Length	Width	Height
400 mm	550 mm	400 mm

VETERINARY EXAMINATION COUCH

MLV-01



Length	Width	Height
1200 mm	600 mm	500-1000 mm

- Adjustable Height using Foot Switch
- Equipped with Waste Drain
- Equipped with Animal Straps



CURA **FLEX**
by **HOSPIcare**

CONTACT US



 +966 5067 05928  www.hospicare.com  info@hospicare.com

 Industrial City District, Medina, Saudi Arabia