

PRIME *Care*

IN-PATIENT CONVENIENCE



TABLE OF CONTENTS

01.	Multi-Functional Prime Care Bed	P. 02
02.	3-Function Regular Side-Rail Bed	P. 04
03.	3-Function CareLuxe Bed	P. 04
04.	3-Function Prime Care Bed	P. 06
05.	3-Function Bariatric Prime Care Bed	P. 08
06.	3-Function Bariatric Prime Care Bed	P. 09
07.	2-Function Manual Bed <i>Regular Side-Rail</i>	P. 10
08.	2-Function Manual Bed <i>With Integrated Toilet</i>	P. 10
09.	Hydraulic Table	P. 11
10.	Over-Bed Table	P. 11
11.	Large Dining Table	P. 11
12.	Bedside Cabinet	P. 11



MULTI-FUNCTIONAL PRIME CARE BED

MLP-01



Length	Width	Height
2000 mm	900 mm	480-810 mm

- Stainless steel I.V Pole
- High Endurance Noiseless Castors; along with Strong Brakes
- Detachable ABS Head and Footboard Supplied with Bumpers
- Optional Accessories:
 - Bed Extension
 - Under Bed Light
 - X-Ray Translucent
 - Urine Bag Holder
 - Battery Backup
 - Central Braking System
 - CPR System

PRIME

 HOSPIcare

 HOSPIcare

3-FUNCTION PRIME CARE BED

MLP-02



Length	Width	Height
2000 mm	900 mm	430-780 mm



- Noiseless Castors with Brakes
- Detachable ABS Head and Footboard
- Stainless Steel I.V. Pole
- Optional Accessories:
 - Battery Backup
 - CPR System
 - X-Ray Translucent
 - Under Bed Light
 - Bed Extension
 - Central Braking System



3-FUNCTION REGULAR SIDE-RAIL BED

MLP-03



Length	Width	Height
2000 mm	900 mm	360-710 mm



- Noiseless Castors with Brakes
- Detachable ABS Head and Footboard
- Optional Accessories:
 - Battery Backup
 - CPR System
 - Over-Bed Table
 - Under Bed Light
 - Bed Extension
 - Central Braking System

3-FUNCTION CARE

MLP-21



- Designed to Prioritize Comfort and Care
- Wooden Headboard and Footboard for a Home-Like Aesthetic
- Noiseless Castors with Brakes for Enhanced Mobility and Safety

- Optional Accessories:
 - Battery Backup
 - CPR System
 - Under Bed Light
 - Bed Extension
 - Central Braking System
 - Over-Bed Table

Length	Width	Height
2000 mm	900 mm	360-710 mm





3-FUNCTION BARIATRIC PRIME CARE

MLP-11



Length	Width	Height
2000 mm	1200 mm	350-610 mm

- Built to Withstand up to **300 KG**
- Noiseless Castors with Brakes
- Detachable ABS Head and Footboard
- Optional Accessories:
 - Battery Backup
 - CPR System
 - X-Ray Translucent
 - Under Bed Light
 - Bed Extension
 - Central Braking System

3-FUNCTION BARIATRIC PRIME CARE BED

MLP-12



Length	Width	Height
2000 mm	1200 mm	350-610 mm

- Built to Withstand up to **300 KG**
- Noiseless Castors with Brakes
- Detachable ABS Head and Footboard
- Stainless Steel I.V Pole
- Optional Accessories:
 - Battery Backup
 - CPR System
 - X-Ray Translucent
 - Under Bed Light
 - Bed Extension
 - Central Braking System

2-FUNCTION MANUAL BEDS

BASIC

MLP-31



Length	Width	Height	Backrest	Legrest
2100 mm	950 mm	530 mm	75°	50°



- I.V Pole
- Urine Bag Holder
- Detachable ABS Head and Footboard

WITH INTEGRATED TOILET

MLP-32



Toilet Base

Length	Width	Height	Backrest	Legrest
2100 mm	950 mm	530 mm	75°	50°

- Coated with electrostatic layer that resists rust and corrosion
- Toilet base movement are carried out by manual crank for convenient use
- I.V pole

TABLES & CABINET

HYDRAULIC

MLH-31



Length	Width	Height
900 mm	450 mm	700-1000 mm

- ABS Top Plate
- Adjustable Height using Gas Spring

OVER-BED

MLH-34



Length	Width
1000 mm	320 mm

- Made of ABS Material
- Light-Weight to Enable Portability

LARGE DINING

MLH-32



Length	Width	Height
900 mm	450 mm	900-1200 mm

- ABS Top Plate
- Height is Manually-Adjusted using an Internal Spring with Knob

BEDSIDE CABINET

MLH-41



Length	Width	Height
430 mm	480 mm	750 mm

- Made of Anti-Bacterial ABS Material
- Lightly Weighed for Ease of Movement



PRIME
Care
by HOSPIcare

CONTACT US



 +966 5067 05928  www.hospicare.com  info@hospicare.com

 Industrial City District, Medina, Saudi Arabia