

PRODUCT DATASHEET


Restiva Titan Pressure-Care Mattress

RESTIVA TITAN · RESTIVA — Pressure Care

Pressure-redistribution mattress combining airflow-enhancing eggshell foam, gel pressure pads and a high-density foam base for pressure ulcer prevention and care.

KEY FEATURES

- Pressure ulcer prevention and treatment
- Enhanced air circulation for patient cooling
- Even pressure distribution
- Honeycomb top layer for better airflow
- Gel pads for pressure redistribution
- High-density foam base for structural support



TECHNICAL SPECIFICATIONS

Available sizes

190 × 90 × 10 cm / 190 × 90 × 14 cm / 200 × 90 × 14 cm