



PRODUCT CATALOG

CONTINUOUS DEVELOPMENT
♣ ULTIMATE SATISFACTION



ICU ELECTRIC BED WITH SIDE CONTROL

MLH-01S

Length	Width	Height	Backrest	Legrest	Trendelenburg	Reverse Trendelenburg
2160 mm	1000 mm	490-730 mm	75°	50°	12°	-12°



- Coated with electrostatic layer that resists rust and corrosion
- Bed movement can be controlled by using remote control and nurse control panel
- Stainless steel I.V Pole • Urine Bag Holder
- Movable on high endurance 125mm noiseless castors; along with strong brakes that stabilize the bed on ground
- Detachable ABS head and footboard supplied with bumpers
- Optional Accessories:
 - Battery Backup - CPR System - Central Braking System - X-Ray Translucent
 - Under Bed Light - Under Bed Tray - Bed Extension

ICU ELECTRIC BED

ABS SIDERAIL



Length	Width	Height	Backrest	Legrest	Trendelenburg	Reverse Trendelenburg
2160 mm	1000 mm	490-730 mm	75°	50°	12°	-12°

- Made of high-endurance materials to ensure safe use
- Coated with electrostatic layer that resists rust and corrosion
- Stainless steel I.V Pole
- Urine Bag Holder
- Movable on high endurance 125mm noiseless castors; along with strong brakes that stabilize the bed on ground
- Detachable ABS head and footboard supplied with bumpers
- Optional Accessories:
 - Bed Extension
 - Under Bed Light
 - X-Ray Translucent
 - Battery Backup
 - Central Braking System
 - CPR System

ICU ELECTRIC BED

REGULAR SIDERAIL

Length	Width	Height	Backrest	Legrest	Trendelenburg	Reverse Trendelenburg
2100 mm	950 mm	470-710 mm	75°	50°	12°	-12°

MLH-01E



- Made with high endurance materials to ensure safe use
- Coated with electrostatic layer that resists rust and corrosion
- Stainless steel I.V Pole
- Urine Bag Holder
- Moves on 125mm diameter noiseless castors with brake
- Aluminum Side-rails
- Detachable ABS Head and Footboard



- Optional Accessories:
 - Battery Backup
 - CPR System
 - Central Braking System
 - Under Bed Light
 - X-Ray Translucent
 - Bed Extension



THREE-FUNCTION ELECTRIC BED WITH SIDE CONTROL

MLH-02S



Length	Width	Height	Backrest	Legrest
2160 mm	1000 mm	490-730 mm	75°	50°

- Coated with electrostatic layer that resists rust and corrosion
- Stainless steel I.V pole
- Moves on 125mm diameter noiseless castors with brakes
- Detachable ABS head and footboard supplied with bumpers
- Optional Accessories:
 - Battery Backup
 - X-ray Translucent
 - Under Bed Light
 - CPR System
 - Central Braking System
 - Bed Extension



THREE-FUNCTION ELECTRIC BED ABS SIDERAIL

MLH-02



- Coated with electrostatic layer that resists rust and corrosion
- Stainless steel I.V pole
- Moves on 125mm diameter noiseless castors with brakes
- Detachable ABS head and footboard supplied with bumpers
- Optional Accessories:
 - Battery Backup
 - Central Braking System
 - Bed Extension
 - CPR System
 - Under Bed Light
 - X-Ray Translucent

Length	Width	Height	Backrest	Legrest
2160 mm	1000 mm	490-730 mm	75°	50°

THREE-FUNCTION ELECTRIC BED REGULAR SIDERAIL

MLH-03



Length	Width	Height	Backrest	Legrest
2100 mm	950 mm	470-710 mm	75°	50°

- Coated with electrostatic layer that resists rust and corrosion
- Stainless Steel I.V pole
- Urine bag holder
- Aluminum side rails
- Moves on 125mm diameter noiseless castors with brake
- Detachable ABS head and footboard
- Optional Accessories:
 - Battery Backup
 - CPR System
 - X-Ray Translucent
 - Under Bed Light
 - Bed Extension
 - Central Braking System



THREE-FUNCTION ELECTRIC BED

REGULAR SIDERAIL | ULTRA LOW HEIGHT

MLH-03E



Length	Width	Height	Backrest	Legrest
2100 mm	950 mm	340-690 mm	75°	50°

- Coated with electrostatic layer that resists rust and corrosion
- Adjustable height reaching very low levels for the ease of use for patients.
- Stainless steel I.V pole
- Urine bag holder
- Aluminum side rails

- Moves on 125mm diameter noiseless castors with brake
- Detachable ABS head and footboard
- Optional Accessories:
 - Battery Backup
 - CPR System
 - Under Bed Light
 - X-Ray Translucent
 - Bed Extension
 - Central Braking System

TWO-FUNCTION ELECTRIC BED ABS SIDERAIL

MLH-04



Length	Width	Height	Backrest	Legrest
2160 mm	1000 mm	550 mm	75°	50°

- Stainless steel I.V pole
- Urine bag holder
- Optional Accessories:
 - X-ray Translucent
 - Backup Battery
 - Central Braking System
 - Under Bed Light
 - Bed Extension

TWO-FUNCTION ELECTRIC BED REGULAR SIDERAIL

MLH-04E



Length	Width	Height	Backrest	Legrest
2100 mm	950 mm	530 mm	75°	50°

- Coated with electrostatic layer that resists rust and corrosion
- Stainless steel I.V pole
- Urine bag holder
- Detachable ABS head and footboard

TWO-FUNCTION ELECTRIC BED WITH CONTROL BUTTONS



Length	Width	Height	Backrest	Legrest
2100 mm	950 mm	530 mm	75°	50°

- Coated with electrostatic layer that resists rust and corrosion
- Bottom-Controlled
- Aluminum side rails
- Stainless steel I.V pole

ONE-FUNCTION ELECTRIC BED WITH CONTROL BUTTON



Length	Width	Height	Backrest
2100 mm	950 mm	530 mm	75°

- Detachable ABS head and footboard
- Urine bag holder

OPTIONAL ACCESSORIES FOR ALL ELECTRIC BEDS

BED EXTENSION



X-RAY TRANSLUCENT BACKREST



UNDER BED TRAY



CENTRAL BRAKING SYSTEM



THREE-FUNCTION MANUAL BED ABS SIDERAIL

MLH-11



Length	Width	Height	Backrest	Legrest
2160 mm	1000 mm	490-730 mm	75°	50°

- Coated with electrostatic layer that resists rust and corrosion
- Stainless steel I.V pole
- Urine bag holder
- Detachable ABS head and footboard supplied with bumpers



THREE-FUNCTION MANUAL BED REGULAR SIDERAIL



Length	Width	Height	Backrest	Legrest
2100 mm	950 mm	470-710 mm	75°	50°



- Coated with electrostatic layer that resists rust and corrosion
- Made with high endurance materials to ensure safe use
- I.V pole
- Urine bag holder
- Detachable ABS head and footboard

TWO-FUNCTION MANUAL BED ABS SIDERAIL



Length	Width	Height	Backrest	Legrest
2160 mm	1000 mm	550 mm	75°	50°



- Coated with electrostatic layer that resists rust and corrosion
- Stainless Steel I.V pole
- Urine bag holder
- Moves on a 125mm diameter noiseless castors with brake
- Detachable ABS head and footboard

TWO-FUNCTION MANUAL BED REGULAR SIDERAIL



Length	Width	Height	Backrest	Legrest
2100 mm	950 mm	530 mm	75°	50°

- Coated with electrostatic layer that resists rust and corrosion
- Made with high endurance materials to ensure safe use
- I.V pole
- Urine bag holder
- Detachable ABS head and footboard

TWO-FUNCTION MANUAL BED

CONVENTIONAL **MLH-12A**



Length	Width	Height	Backrest	Legrest
1980 mm	890 mm	480 mm	75°	50°

- An economic bed
- Made with high-endurance materials to ensure safe use
- Coated with electrostatic layer that resists rust and corrosion
- Manually controlled bed using a crank
- I.V pole

MLH-12B WITH INTEGRATED TOILET



Length	Width	Height	Backrest	Legrest
2100 mm	950 mm	530 mm	75°	50°

- Coated with electrostatic layer that resists rust and corrosion
- Toilet base movement are carried out by manual crank for convenient use
- I.V pole

ONE-FUNCTION MANUAL BED

REGULAR SIDERAIL

MLH-13



Length	Width	Height	Backrest
2100 mm	950 mm	530 mm	75°

- Coated with electrostatic layer that resists rust and corrosion
- Made with high-endurance materials to ensure safe use
- I.V pole
- Detachable ABS head and footboard

MLH-13A

CONVENTIONAL



Length	Width	Height	Backrest
1980 mm	890 mm	480 mm	75°

- An economic bed
- Coated with electrostatic layer that resists rust and corrosion
- Provided with an I.V pole



THREE-FUNCTION ELECTRIC HEMODIALYSIS CHAIR

Length	Width	Height
2050 mm	600 mm	600-750 mm

MLH-06



- Coated with electrostatic layer that resists rust and corrosion
- Can be extended into a full bed
- Movable armrests
- Adjustable footrest
- Adjustable height
- Provided with lower shelf and paper roll holder
- Provided with an adjustable stainless steel I.V pole
- Moves on a 125mm diameter noiseless castors with brake
- Optional Accessories:
 - Over bed table

TWO-FUNCTION ELECTRIC HEMODIALYSIS CHAIR

TWO-FUNCTION MANUAL HEMODIALYSIS CHAIR

Length	Width	Height
2050 mm	600 mm	600 mm

MLH-07



- Coated with electrostatic layer that resists rust and corrosion
- Can be extended into a full bed
- Movable armrests
- Adjustable footrest
- Provided with lower shelf and paper roll holder



MLH-15



- Provided with an adjustable stainless-steel I.V pole
- Moves on 125mm diameter noiseless castors with brakes
- Optional Accessories:
 - Over bed table

ONE-FUNCTION PEDIATRIC BED

Length	Width	Height	Backrest
1270 mm	600 mm	1200 mm	75°

REGULAR

MLH-21



FIXED

MLH-21E



- Safe for children aging 2-8 years
- Easy elevating and lowering side rails and easy tilting backrest

INFANT BED

WITH DRAWER

MLH-22



Length	Width	Height
820 mm	460 mm	960 mm

- Painted with electrostatic layer that resists rust and corrosion
- Moves on a 75mm diameter castors with brakes
- Manual trendelenburg and reverse
- The bed has drawer and basket for extra storage of baby items

MLH-23



Length	Width	Height
820 mm	460 mm	920 mm

- Manual one step trendelenburg
- Moves on a 75mm diameter castors with brakes

OVER BED TABLES

HYDRAULIC

MLH-31



- ABS top plate
- Adjustable height using gas spring

Length	Width	Height
900 mm	450 mm	700-1000 mm

MLH-34

OVER-BED



Length	Width
880 mm	320 mm

- The table is made of ABS material
- Light-weight table to enable portability
- Adjustable to the bed width

LARGE DINING

MLH-32



- ABS top plate
- Coated with electrostatic technology and Nickel-Chrome to resist rust and corrosion
- Height is manually-adjusted using an internal spring with knob
- Moves on a 75mm diameter castors with brakes

Length	Width	Height
900 mm	450 mm	900-1200 mm

SMALL DINING

MLH-33



- ABS top plate
- Coated with electrostatic technology and Nickel-Chrome to resist rust and corrosion
- Height is manually-adjusted
- Moves on a 50mm diameter castors with brakes

Length	Width	Height
810 mm	320 mm	820-1150 mm

BEDSIDE CABINETS

ABS

MLH-41



Length	Width	Height
430 mm	480 mm	750 mm

- Made of anti-bacterial ABS material
- Moveable with the four, high quality castors
- Lightly weighed and can be moved easily

MLH-46

STAINLESS STEEL



Length	Width	Height
430 mm	330 mm	630 mm

- Stainless steel top plate
- Coated with electrostatic technology to resist rust and corrosion
- Based on 4 adjustable bases for more stability
- Light- weighted

I.V. INFUSION STAND

MLH-42



Height
1100-2100 mm

- Made of stainless steel
- 4 Hooks
- Height is manually-adjusted
- Moves on plastic base with five 50mm diameter

BED MATTRESSES

- Made of high density sponge to achieve optimum comfort for the patient
- Coated by leather to prevent leakage of liquids and for ease of cleaning



4-SECTION MATTRESS **MLH-43**



Length	Width	Height
1980 mm	890 mm	100 mm

2-SECTION MATTRESS **MLH-44**



Length	Width	Height
1980 mm	890 mm	100 mm

1-PIECE MATTRESS

MLH-45

Length	Width	Height
1980 mm	890 mm	100 mm

MATTRESS COLORS

(OTHER COLORS ARE AVAILBLE UPON REQUEST)



APPAREL'S CABINET

- Made of electrostatic-coated steel
- Locker is mounted on four rubber footings

1-DOOR SIDE

MLH-51



Length	Width	Height
450 mm	400 mm	1800 mm

2-DOOR SIDE

MLH-52



Length	Width	Height
900 mm	400 mm	1800 mm

LOCKER | 3-DOOR SIDE

MLH-53



Length	Width	Height
450 mm	400 mm	1800 mm

EXAMINATION COUCH

MANUAL

MLC-01




Length	Width	Height
1850 mm	600 mm	820 mm

- Coated with electrostatic layer and nickel-chrome to resist rust and corrosion, and easily sanitized to guarantee longer life span
- Manual backrest tilting
- Accessories:
 - Paper roll holder

EXAMINATION COUCH

FIXED

MLC-01E



Length	Width	Height
1850 mm	600 mm	820 mm

- Coated with electrostatic layer and nickel-chrome to resist rust and corrosion, and easily sanitized to guarantee longer life span
- Accessories:
 - Paper roll holder

GYNAECOLOGY BED

MLC-02



Length	Width	Height
1850 mm	600 mm	820 mm

- Coated with electrostatic layer and nickel-chrome to resist rust and corrosion, and easily sanitized to guarantee longer life span
- Manual back and knee rest tilting
- Can be extended into a full bed
- Equipped with stainless steel basin and adjustable leg holders
- Based on 4 plastic plugs

MULTI-FUNCTIONAL CHAIR

MLC-03



Length	Width	Height
1750 mm	550 mm	470 mm

- Coated with electrostatic layer that resists rust and corrosion
- Can be extended into a full bed
- Movable Stainless Steel armrests
- *Optional Accessories:*
 - I.V Infusion Stand
 - Arm Holder

ELECTRIC EXAMINATION COUCH

MLC-05



- Coated with electrostatic layer that resists rust and corrosion
- Adjustable height using foot switch
- Equipped with a paper roll holder
- Based on 4 adjustable bases for more stability

Length	Width	Height	Backrest Tilting
2000 mm	600 mm	530-1000 mm	75°

ELECTRIC DELIVERY BED

MLC-06



Length	Width	Height
1850 mm	600 mm	630-980 mm

- Coated with electrostatic layer that resists rust and corrosion
- Adjustable back and kneerest tilting
- The movement of the back and foot is controlled by a multi-height serrated ruler
- Adjustable height using foot switch

- Can be extended into a full bed
- Equipped with a stainless-steel basin and adjustable leg holders
- Based on 4 adjustable bases for more stability
- Optional Accessory: Siderails

MLC-11

DESK

Length	Width	Height
1200 mm	600 mm	800 mm

- Provided with a single drawer
- Top of the desk is made of high quality MDF
- Mounted on rubber footings for safety and higher stability

MLC-12

RECEPTION COUNTER

- Lower Plate Dimensions: 1200mm * 600mm
- Upper Plate Dimensions: 900mm * 250mm
- Provided with a single drawer
- Pedals are added in the lower part of the counter to offer the highest level of comfort

3-DESK DRAWERS

MLC-13

Length	Width	Height
400 mm	400 mm	650 mm

- Mounted on castors
- Optional accessory for item MLC-11 and MLC-12

WAITING ROOM SEATS

AVAILABLE IN 2, 3 & 4 UNIT SEATS



- Made of electrostatic-painted steel with high durability
- Seat made of high-density foam to provide optimum comfort



- Ergonomically designed to provide optimum comfort
- Plastic seats with smooth texture for optimum comfort and ease of cleaning

WAITING ROOM CHAIR



- Made of electrostatic painted steel with anti-rusting property
- Seat made of high density foam with leather cover to provide optimum comfort

COLORS (OTHERS AVIALBLE UPON REQUEST)



EXAMINATION CHAIRS

- Hydraulic handle under the chair's seat controls elevating and lowering of the chair
- 360° rotatable seat
- High-density sponge leather seat base to provide physician's optimum comfort



INSTRUMENTS' TABLE

MLC-31



Length	Width	Height
400 mm	600 mm	800 mm

- Stainless-steel structure
- Moves on a 50mm diameter castors

MAYO INSTRUMENTS' TABLE

MLC-34



Length	Width	Height
460 mm	630 mm	830-1170 mm

- Steel structure with chrome plating to prevent rust and corrosion
- Stainless-steel removable tray
- Adjustable height using a hand knob
- Moves on a 75mm diameter castors with brakes

X-RAY FILM VIEWER

1 SCREEN

Length	Width
400 mm	450 mm

MLC-41

3 SCREENS

Length	Width
1200 mm	450 mm

MLC-43



2 SCREENS

Length	Width
800 mm	450 mm

MLC-42

4 SCREENS

Length	Width
1600 mm	450 mm

MLC-44

- Body made of electrostatic-coated steel and front surface high quality light diffuser

INSTRUMENTS' CABINETS

2-SLIDING DOORS

MLC-51



Length	Width	Height
900 mm	350 mm	1600 mm

- Made of steel, plated with electrostatic material to prevent rust and corrosion
- Two sliding doors that makes it easy to use
- Doors have transparent glass
- The lower part of the wardrobe has 2 doors for storage, accompanied by a space for applying a lock for safety

ONE DOOR

MLC-52



Length	Width	Height
450 mm	350 mm	1800 mm

TWO DOORS

MLC-53



Length	Width	Height
900 mm	350 mm	1800 mm

WARD SCREENS

3-SECTION

MLC-61



Length	Width
1800 mm	1500 mm

4-SECTION

MLC-61Q

- Chrome-plated steel tubular frame to prevent rust and corrosion
- Mounted on 50mm castors for easy movement
- Also available in distinguished fabric with code **MLC-61L**

BOWL STANDS

MLC-62



MLC-63



- Made of chrome-plated steel that prevents rusting and corrosion
 - Bowls can be easily removed and installed
- Five spoke heavy duty plastic base with castors

STEPS

- Covered with leather that avoids sliding

MLC-64 SINGLE STEP



Length	Width	Height
400 mm	400 mm	250 mm

DOUBLE STEP **MLC-65**



Length	Width	Height
400 mm	550 mm	400 mm

1-FUNCTION EMERGENCY TROLLEY

Length	Width	Height	Backrest Tilting	OR	Length	Width	Height	Backrest Tilting
1750 mm	600 mm	700 mm	75°		1900 mm	600 mm	700 mm	75°



MLT-02

- Back section's angle is adjusted through hand crank
- Coated with electrostatic layer that resists rust and corrosion
- Fixed height
- Provided with I.V pole and cylinder holder
- Aluminum side rail
- Moves on a 125mm diameter noiseless castors with brake



MLT-03



2-FUNCTION EMERGENCY TROLLEY

- Coated with electrostatic layer that resists rust and corrosion
- Adjustable height from 500 to 700 mm
- Provided with I.V pole and cylinder holder
- Aluminum side rail
- Moves on a 125mm diameter noiseless castors with brakes

STRETCHER TROLLEY

WITH SEPARABLE STRAPS



MLT-04

Length	Width	Height
2000 mm	630 mm	780 mm

- Top plate is made of stainless steel
- Easily detachable stretcher
- Equipped with straps for patient safety
- Coated with electrostatic technology that resists rust and corrosion
- Moves on a 125mm diameter castors with brake

PATIENT TROLLEY

SEPARABLE



MLT-05


Length	Width	Height
1950 mm	530 mm	800 mm

- The stretcher can be detached from the trolley and re-attached easily
- The stretcher is padded with high density foam and covered with leather, for optimum comfort for the patient
- Manual backrest tilting
- Stretcher has legs to stabilize it on the floor

CRASH TROLLEY

Length	Width	Height
840 mm	500 mm	950 mm

- Consists of five drawers, divided from inside
- Accessories illustrated on the next diagram



MLT-11



STAINLESS-STEEL CRASH TROLLEY **MLT-12**



Length	Width	Height
600 mm	500 mm	900 mm

- Stainless-steel structure
- Moves on a 100 mm diameter castors with brakes
- Provided with:
 - Adjustable I.V. Pole
 - Oxygen cylinder holder
 - ABS cardiac boards
 - Defibrillator Shelf
 - Waste bin



FOLDABLE LAUNDRY TROLLEY



Length	Width	Height
600 mm	600 mm	1000 mm

MLT-13

- Equipped with a detachable sack for collecting dirty linen
- Movable on 4 castors
- Equipped with 2 hooks on the side for holding bags
- Can be easily folded and stored
- Fixed laundry trolley is also available with code **MLT-13O**

ROOM SERVICE TROLLEY



Length	Width	Height
1000 mm	500 mm	1000 mm

MLT-14


- Equipped with 3 shelves
- Equipped with a detachable sack for collecting dirty linen
- Movable on 4 castors



CONTACT US:

HQ: 14-12 Rond Point des Champs-Élysées, 75008, Paris, France
 Factory KSA: Industrial City District, Medina, Saudi Arabia
Factory Egypt: 171,173,175 First Industrial Zone, Badr City, Cairo, Egypt

 m.swilam@hospicare.fr
 www.hospicare.fr

 +337 84041864
+966 50 670 5928
+202 27933393/4/6