

# CURA <sup>^</sup>FLEX

OUT-PATIENT SOLUTIONS



# TABLE OF CONTENTS

<b>01.</b> Electric Examination Couch	<b>P. 03</b>
<b>02.</b> 2-Function Electric Examination Couch	<b>P. 04</b>
<b>03.</b> 3-Section Electric Examination Couch	<b>P. 05</b>
<b>04.</b> Fixed Examination Couch	<b>P. 06</b>
<b>05.</b> Manual Examination Couch	<b>P. 06</b>
<b>06.</b> Electric Gynaecology Examination Couch	<b>P. 07</b>
<b>07.</b> Manual Gynaecology Couch	<b>P. 08</b>
<b>08.</b> 4-Function Hemodyalisis Chair	<b>P. 10</b>
<b>09.</b> 3-Function Trendelenburg Hemodyalisis Chair	<b>P. 11</b>
<b>10.</b> 3-Function Hemodyalisis Chair	<b>P. 12</b>
<b>11.</b> Manual 2-Function Hemodyalisis Chair	<b>P. 13</b>
<b>12.</b> Electric 2-Function Hemodyalisis Chair	<b>P. 13</b>
<b>13.</b> Multi-Functional Chair	<b>P. 14</b>
<b>14.</b> Instruments' Tables	<b>P. 15</b>
<b>15.</b> Ward Screens   Bowl Stands   Steps	<b>P. 16</b>
<b>16.</b> Veterinary Examination Couch	<b>P. 17</b>



# ELECTRIC EXAMINATION COUCH

MLC-05



Foot Switch

Length	Width	Height	Backrest Tilting
1850 mm	600/650 mm	570-900 mm	80°

- Equipped with a Paper Roll Holder
- Height Adjustable by Electric Foot Switch
- Optional Accessories:
  - I.V. Stand
  - Headrest Opening

## 2-FUNCTION ELECTRIC EXAMINATION COUCH

MLC-05S



Length	Width	Height	Backrest Tilting
1850 mm	600/650 mm	570-900 mm	80°

- Equipped with a Paper Roll Holder
- Backrest and Height Adjustable by Electric Hand Control
- Optional Accessories:
  - I.V. Stand
  - Headrest Opening

## 3-SECTION ELECTRIC EXAMINATION COUCH

MLC-04



Foot Switch

Length	Width	Height	Backrest Tilting
1750 mm	600/650 mm	600-950 mm	80°

- Equipped with a Paper Roll Holder
- Height Adjustable by Electric Foot Switch
- Optional Accessories:
  - I.V. Stand
  - Headrest Opening

# EXAMINATION COUCHES

## FIXED EXAMINATION COUCH

**MLC-01E**



Length	Width	Height
1850 mm	600 mm	820 mm

- Easily sanitized to guarantee longer life span
- Accessories:
  - Paper Roll Holder

## MANUAL EXAMINATION COUCH

**MLC-01**



Length	Width	Height
1850 mm	600/650 mm	700/820 mm

- Easily sanitized to guarantee longer life span
- Manual Backrest Tilting
- Accessories:
  - Paper Roll Holder
  - I.V. Pole

# ELECTRIC GYNAECOLOGY EXAMINATION COUCH

MLC-06



Adjustable Leg Holders



S.S Basin

Foot Switch



Length	Width	Height
1750 mm	600/650 mm	600-950 mm

- Adjustable back and knee rest tilting
- Adjustable height using foot switch
- Can be extended into a full bed
- Optional Accessory:
  - Side Rails
  - I.V. Pole



# MANUAL GYNAECOLOGY COUCH

MLC-02



Length	Width	Height
1750 mm	600/650 mm	700/820 mm

- Easily Sanitized to Guarantee Longer Life Span
- Manual Back and Knee Rest Tilting
- Can be Extended into a Full Bed
- Optional Accessory:
  - I.V. Pole
  - Paper Roll Holder



# 4-FUNCTION HEMODYALYSIS CHAIR

MLH-061

Padded  
Seat Cushion

Over-Seat  
Table



Length	Width	Height	Trendelenburg
1950 mm	600 mm	600-800 mm	16°



- Can be Extended into a Full Bed
- Movable Armrests
- Adjustable Footrest
- Provided with Lower Shelf
- **Optional Accessories:**
  - Paper Roll Holder
  - Over-Seat Table

## 3-FUNCTION TRENDELENBURG HEMODYALYSIS CHAIR

MLH-06T



Length	Width	Height	Trendelenburg
1950 mm	600 mm	550 mm	16°



- Can be Extended into a Full Bed
- Emergency Trendelenburg Function
- Movable Armrests
- Adjustable Footrest
- Provided with Lower Shelf
- **Optional Accessories:**
  - Paper Roll Holder
  - Over-Seat Table

# 3-FUNCTION HEMODIALYSIS CHAIR

MLH-06



Length	Width	Height
1950 mm	600 mm	550-750 mm

- Can be Extended into a Full Bed
- Movable Armrests
- Adjustable Footrest
- Adjustable Height
- Provided with Lower Shelf
- Optional Accessories:
  - Over Bed Table
  - Paper Roll Holder

# HEMODIALYSIS CHAIRS

## MANUAL 2-FUNCTION HD CHAIR

## MLH-15



Length	Width	Height
1950 mm	600 mm	550 mm

- Provided with an Adjustable Stainless Steel I.V Pole
- Noiseless Castors with Brakes
- Optional Accessories:
  - Over Seat Table

## ELECTRIC 2-FUNCTION HD CHAIR

## MLH-07



- Can be Extended into a Full Bed
- Movable Armrests
- Adjustable Footrest
- Provided with Lower Shelf
- Optional Accessories:
  - Paper Roll Holder
  - Over Seat Table

Length	Width	Height
1950 mm	600 mm	550 mm

# MULTI-FUNCTIONAL CHAIR

MLC-03



Length	Width	Height
1750 mm	550 mm	470 mm

- Can be Extended into a Full Bed
- Movable Stainless Steel Armrests
- Optional Accessories:
  - I.V Infusion Stand
  - Arm Holder

# INSTRUMENTS' TABLES

**DUAL-TIER TABLE**

**MLC-31**



Length	Width	Height
400/500 mm	600/700 mm	800 mm

- Stainless Steel Structure
- Moves on a 50mm Diameter Castors

**TRIPLE-TIER TABLE**

**MLC-32**



Length	Width	Height
400/500 mm	600/700 mm	800 mm

- Stainless Steel Structure
- Moves on a 50mm Diameter Castors

**MAYO TABLE**

**MLC-34**



Length	Width	Height
460 mm	630 mm	830-1170 mm

- Steel Structure with Chrome Plating
- Stainless Steel Removable Tray
- Adjustable Height using a Hand Knob
- 75mm Diameter Castors with Brakes

**WITH DRAWERS**

**MLC-31D**



Length	Width	Height
400/500 mm	600/700 mm	800 mm

- Stainless Steel Structure
- 2 Stainless Steel Drawers
- Moves on a 50mm Diameter Castors



**WARD SCREENS**    **MLC-61**



Length	Width
1800 mm	1500 mm

- Chrome-Plated Steel Tubular Frame
- Available in Distiguated Fabric (MLC-61L)
- Available in 4-Section (MLC-61Q)

**MLC-62**    **BOWL STANDS**    **MLC-63**



- Chrome-Plated Steel
- Easily Removed and Installed
- Heavy Duty Base with Castors



**MLC-64**    **STEPS**    **MLC-65**

**SINGLE**



Leather-Covered  
Anti-Slip Protection

**DOUBLE**



Length	Width	Height
400 mm	400 mm	250 mm

Length	Width	Height
400 mm	550 mm	400 mm

# VETERINARY PROCEDURE TABLE

MLV-01



Length	Width	Height
1200 mm	600 mm	500-1000 mm

- Adjustable Height using Foot Switch
- Equipped with Waste Drain
- Equipped with Animal Straps



**CURA** <sup>FLEX</sup>  
by HOSPIcare

## CONTACT US



 +966 5067 05928  [www.hospicare.com](http://www.hospicare.com)  [info@hospicare.com](mailto:info@hospicare.com)

 Industrial City District, Medina, Saudi Arabia

# PRIME *Care*

**IN-PATIENT CONVENIENCE**



# TABLE OF CONTENTS

<b>01.</b>	Multi-Functional Prime Care Bed	<b>P. 02</b>
<b>02.</b>	3-Function Regular Side-Rail Bed	<b>P. 04</b>
<b>03.</b>	3-Function CareLuxe Bed	<b>P. 04</b>
<b>04.</b>	3-Function Prime Care Bed	<b>P. 06</b>
<b>05.</b>	3-Function Bariatric Prime Care Bed	<b>P. 08</b>
<b>06.</b>	3-Function Bariatric Prime Care Bed	<b>P. 09</b>
<b>07.</b>	2-Function Manual Bed   <i>Regular Side-Rail</i>	<b>P. 10</b>
<b>08.</b>	2-Function Manual Bed   <i>With Integrated Toilet</i>	<b>P. 10</b>
<b>09.</b>	Hydraulic Table	<b>P. 11</b>
<b>10.</b>	Over-Bed Table	<b>P. 11</b>
<b>11.</b>	Large Dining Table	<b>P. 11</b>
<b>12.</b>	Bedside Cabinet	<b>P. 11</b>



## MULTI-FUNCTIONAL PRIME CARE BED

### MLP-01



Length	Width	Height
2000 mm	900 mm	480-810 mm

- Stainless steel I.V Pole
- High Endurance Noiseless Castors; along with Strong Brakes
- Detachable ABS Head and Footboard Supplied with Bumpers
- Optional Accessories:
  - Bed Extension
  - Under Bed Light
  - X-Ray Translucent
  - Urine Bag Holder
  - Battery Backup
  - Central Braking System
  - CPR System

PRIME

 HOSPIcare

 HOSPIcare



3-FUNCTION  
**PRIME CARE  
 BED**  
 MLP-02



Length	Width	Height
2000 mm	900 mm	430-780 mm



- Noiseless Castors with Brakes
- Detachable ABS Head and Footboard
- Stainless Steel I.V Pole
- Optional Accessories:
  - Battery Backup
  - CPR System
  - X-Ray Translucent
  - Under Bed Light
  - Bed Extension
  - Central Braking System



## 3-FUNCTION REGULAR SIDE-RAIL BED

MLP-03



Length	Width	Height
2000 mm	900 mm	360-710 mm



- Noiseless Castors with Brakes
- Detachable ABS Head and Footboard
- Optional Accessories:
  - Battery Backup
  - CPR System
  - Over-Bed Table
  - Under Bed Light
  - Bed Extension
  - Central Braking System

## 3-FUNCTION CARE

MLP-21



- Designed to Prioritize Comfort and Care
- Wooden Headboard and Footboard for a Home-Like Aesthetic
- Noiseless Castors with Brakes for Enhanced Mobility and Safety

- Optional Accessories:
  - Battery Backup
  - CPR System
  - Under Bed Light
  - Bed Extension
  - Central Braking System
  - Over-Bed Table


Length	Width	Height
2000 mm	900 mm	360-710 mm





## 3-FUNCTION BARIATRIC PRIME CARE

MLP-11



Length	Width	Height
2000 mm	1200 mm	350-610 mm

- Built to Withstand up to **300 KG**
- Noiseless Castors with Brakes
- Detachable ABS Head and Footboard
- Optional Accessories:
  - Battery Backup
  - CPR System
  - X-Ray Translucent
  - Under Bed Light
  - Bed Extension
  - Central Braking System

3-FUNCTION  
BARIATRIC  
PRIME CARE BED

MLP-12



Length	Width	Height
2000 mm	1200 mm	350-610 mm

- Built to Withstand up to **300 KG**
- Noiseless Castors with Brakes
- Detachable ABS Head and Footboard
- Stainless Steel I.V Pole
- Optional Accessories:
  - Battery Backup
  - CPR System
  - X-Ray Translucent
  - Under Bed Light
  - Bed Extension
  - Central Braking System

## 2-FUNCTION MANUAL BEDS

**BASIC** **MLP-31**



Length	Width	Height	Backrest	Legrest
2100 mm	950 mm	530 mm	75°	50°

- I.V Pole
- Urine Bag Holder
- Detachable ABS Head and Footboard

**WITH INTEGRATED TOILET** **MLP-32**



Length	Width	Height	Backrest	Legrest
2100 mm	950 mm	530 mm	75°	50°

- Coated with electrostatic layer that resists rust and corrosion
- Toilet base movement are carried out by manual crank for convenient use
- I.V pole

# TABLES & CABINET

HYDRAULIC

MLH-31



Length	Width	Height
900 mm	450 mm	700-1000 mm

- ABS Top Plate
- Adjustable Height using Gas Spring

OVER-BED

MLH-34



Length	Width
1000 mm	320 mm

- Made of ABS Material
- Light-Weight to Enable Portability

LARGE DINING

MLH-32



Length	Width	Height
900 mm	450 mm	900-1200 mm

- ABS Top Plate
- Height is Manually-Adjusted using an Internal Spring with Knob

BEDSIDE CABINET

MLH-41



Length	Width	Height
430 mm	480 mm	750 mm

- Made of Anti-Bacterial ABS Material
- Lightly Weighed for Ease of Movement





**PRIME**  
*Care*  
by HOSPICARE

## CONTACT US



 +966 5067 05928  [www.hospicare.com](http://www.hospicare.com)  [info@hospicare.com](mailto:info@hospicare.com)

 Industrial City District, Medina, Saudi Arabia